



MRI patient checklist

The system delivering HF10™ is approved for MRI scanning of your head and extremities under certain conditions.

When a doctor requests an MRI scan:

1. Tell the doctor who requested the MRI scan that you have an implanted spinal cord stimulation device.
2. Ask the doctor who implanted your system:
 - Can my system safely undergo the ordered MRI scan?
3. When scheduling your MRI appointment, provide the following information:
 - You have a Nevro spinal cord stimulator.
 - The contact information of the doctor who implanted your device.
 - The MRI guidelines, which can be found at HF10.com/resources/existing-patients.
4. Take your patient ID card to the MRI appointment.
5. Prior to your MRI scan, make sure your stimulation is turned OFF. You can turn stimulation OFF with your Remote Control.
6. After the MRI scan, turn your system back ON with your Remote Control and contact your doctor to ensure that your HF10 is programmed to pre-MRI settings.

Note:

- The system is not currently approved for MRI scans when used with a surgical lead (also known as a "paddle lead") or leads from other manufacturers.
- It's possible that additional scanning options will be available in the future, so please refer back to the HF10.com website to confirm the most current guidelines.

