HFX App Compatibility List

Below are the Apple iPhones and supported OS version that are compatible with the HFX App:

Device Name

iPhone 7	iPhon	iPhone 11 Pro		iPhone 13 Pro
iPhone 7 Plus	iPhon	iPhone 11 Pro Max		iPhone 13 Pro Max
iPhone 8	iPhon	iPhone SE 2nd Gen		iPhone SE 3rd Gen
iPhone 8 Plus	iPhon	iPhone 12 Mini		iPhone 14
iPhone X	iPhon	iPhone 12		iPhone 14 Pro
iPhone XS	iPhon	iPhone 12 Pro		iPhone 14 Plus
iPhone XS Max	iPhon	iPhone 12 Pro Max		iPhone 14 Pro Max
iPhone XR	iPhon	iPhone 13 Mini*		
iPhone 11	iPhon	iPhone 13		
Minimum iOS Version Minimum Available Memory (MB)				
13		5		

Additional Information

57

- If your Apple device isn't listed on this page, it's not supported by the HFX App.
- The HFX App is not compatible with iPad or iPod devices.
- * If you're using an iPhone 13 Mini, please keep your phone within 2 meters (6 ½ feet) from your HFX iQ device when using the HFX App. This will ensure the connection between the two devices is strong enough so your system works properly.

SENZA, SENZA II, SENZA OMNIA, OMNIA, HF10, the HF10 logo, 10 kHz Therapy, HFX, the HFX logo, HFX iQ, the HFX iQ logo, HFX Algorithm, HFX CONNECT, the HFX Connect logo, HFX ACCESS, the HFX Access logo, HFX COACH, the HFX Coach logo, HFX CLOUD, the HFX Cloud logo, RELIEF MULTIPLIED, the X logo, NEVRO, and the NEVRO logo are trademarks or registered trademarks of Nevro Corp. ©2022 Nevro Corp. All Rights Reserved.

ML-00031 Rev. D



HFX_i